# CardioSecur Pro

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The screenshot shows the 'ISP Settings' screen. At the top, it says 'Patients' and 'ISP'. Below the header is a diagram showing a patient with an ECG lead on their chest, connected to a central 'ISP' logo, which is then connected to a doctor sitting at a computer. Below the diagram are three input fields: 'email address' with the value 'doc@cardiosecur.com', 'Password' with masked characters, and 'Synchronisation' with a toggle switch turned on. Below these fields is a section titled 'ISP Registration' with the following text: 'The synchronization with the ISP offers the following options: - You receive all features for CardioSecur PRD that have been assigned to your account. If you wish to add one of the listed features to your subscription, please contact our customer support team via email at info@cardiosecur.com. - Your CardioSecur PRD app serves as a cockpit for all ECG data that has been generated with your patients' CardioSecur ACTIVE devices. You can view, analyze and export this data. Please enter your email address and password and select synchronization. The synchronization will start automatically.' At the bottom, there are two options: '15-lead ECG' and '18-lead ECG', both with checkmarks.

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# CardioSecur Pro

User Manual (EN)  
(with detailed App description)

Version 10.0  
26.10.2017

# Index

1	Symbols	3
2	Intended Use	4
3	Safety Information	5
4	Getting Started	7
5	Data Backup	7
6	ISP	7
7	GDT-Interface	8
	7.1 Overview	8
	7.2 WebDAV-Server	8
	7.3 Data structures	10
	7.4 Process	11
8	Care and Maintenance	12
9	Technical Data and Operating Conditions	13
10	Electromagnetic Compliance (EMC)	14
11	Disposal	16
12	Labelling	16
13	Contact	16
14	iPad	17
	14.1 ECG Recording	17
	14.2 Emergency ECG Recording	19
	14.3 Menu	20
	14.4 Patient Administration	23
	14.5 Patient Details	24
	14.6 ECG Display and Functions	25
	14.7 Archive	28
	14.8 Export	29
15	iPhone	30
	15.1 ECG Recording	30
	15.2 Quick reading	32
	15.3 Menu	33
	15.4 Patient Administration	37
	15.5 Patient Details	38
	15.6 ECG Display and Functions	39
	15.7 Archive	43
	15.8 Export	44

This user manual describes the use of CardioSecur Pro in conjunction with the following iOS devices (iPhone/iPad is a trademark of Apple Inc.):


- iPad 2, iPad (3<sup>rd</sup> generation), iPad 4, iPad Air and iPad Mini.
- iPhone 4S, iPhone 5/5S/5C, iPhone 6/6 plus, iPhone 7/7 plus.

For simplicity herein referred to collectively as iOS device and CardioSecur.

Additional information can be found at [www.cardiosecur.com/en/pro-manual](http://www.cardiosecur.com/en/pro-manual).

# 1 Symbols

Packaging or instructions included in the scope of delivery may display the following symbols (color and size can vary):

Symbol	Meaning
	Use disposable electrodes before the displayed date of expiry.
<b>LOT</b>	Lot denomination
	This product complies with the requirements of the directive for medical devices 93/42/EWG.
<b>REF</b>	Item number
	Please read the user manual
	Warning Warning concerns safety-relevant information
	Information Information concerns details that are worth knowing
	Manufacturer's data
	Store at indicated temperature with unbroken seal.
	For single use only
	Latex-free
	CardioSecur and its accessories - like any other electronic device - contain metal and plastic parts, which need to be disposed of separately after the operating life to prevent a negative environmental impact. Send CardioSecur and its associated accessories with sufficient postage, if possible in its original packaging, to Personal MedSystems GmbH for disposal free of charge.
	Application unit of type BF. An application unit of type BF grants protection against electric shock, in particular with regard to permissible discharge current, according to the requirements of norm EN 60601-1.


## 2 Intended Use

CardioSecur is intended to record, evaluate, document and diagnose a resting electrocardiogram (ECG) for adult males and females.

CardioSecur is a medical electrical system comprising an ECG cable, an iOS device and an app for an iOS device. CardioSecur records a 22-lead ECG using 4 electrodes.

Only medically trained personnel may use CardioSecur.

To record an ECG with CardioSecur, the disposable electrodes need to be positioned on the patient's thorax according to the following diagram:

- White electrode:** At the upper end of the sternum (on the bone)
- Yellow electrode:** At the lower end of the sternum (on the bone)
- Red electrode:** In the middle of the right axillae at the height of the lower end of the sternum (the same height as the yellow electrode)
- Green electrode:** In the middle of the left axillae at the height of the lower end of the sternum (the same height as the yellow electrode)

### 3 Safety Information

Take account of the following warnings and information. In case of your non-compliance you run the risk of causing malfunction, damage or bodily injury.

**⚠ Warning:** CardioSecur is provided for the exclusive use by qualified physicians or personnel under their direct supervision. The numerical and graphical results from a recording must be examined with respect to the patient's overall clinical condition. It is the physicians responsibility to make the diagnosis or to obtain expert opinion on the results and to institute correct treatment if indicated.

**⚠ Warning:** In order to ensure sufficient signal quality and to prevent distorted results, only use disposable electrodes that are unused, have not expired and are at room temperature.

**⚠ Warning:** Remove all disposable electrodes from the skin after a reading, as prolonged adhesion may cause skin irritation.

**⚠ Warning:** Do not use CardioSecur in case of sensitive skin or allergies. This could cause rash, irritation of the skin or eczema.

**⚠ Warning:** Do not use CardioSecur when using an external defibrillator.

**⚠ Warning:** Do not perform readings in places exposing CardioSecur to high electromagnetic radiation (e.g. TV sets, monitors, dish washers).

**⚠ Warning:** Do not expose CardioSecur to any electrostatics. Discharge any electrostatics from your and your patient's body before touching CardioSecur.

**⚠ Warning:** Do not expose CardioSecur to any physical shocks or vibrations. Do not drop it on the floor or step on it.

**⚠ Warning:** Do not place anything onto CardioSecur as this may damage it.

**⚠ Warning:** Do not open, take apart or repair CardioSecur. Do not transform CardioSecur. CardioSecur does not contain any parts that need to be maintained by its user.

**⚠ Warning:** Store and use CardioSecur only under the conditions described in section **Error! Reference source not found., Error! Reference source not found..**

**⚠ Warning:** Do not sterilize CardioSecur in a steam sterilizer (autoclave) or in a gas sterilizer (ethylene oxide, formaldehyde, ozone etc.).

**⚠ Warning:** Do not use CardioSecur in the vicinity of combustible anaesthetics, drugs or pressurized oxygen.

**⚠ Warning:** Only use CardioSecur if it has no apparent damage.

**⚠ Warning:** Note the information in section **Error! Reference source not found., Error! Reference source not found..**

⚠ **Warning:** CardioSecur utilizes an ECG method computing up to 22-lead ECG from 3 leads read with 4 electrodes. Hence, deviations might arise in comparison to a 12-lead ECG read with 10 electrodes.

⚠ **Warning:** Do not perform ECG readings if the skin is wet at the attachment points of the disposable electrodes, e.g. due to sweat or after bathing.

⚠ **Warning:** Proper functioning of the disposable electrodes is influenced by large amounts of body hair, cream or oily skin. Regularly remove hair from your skin at the attachment points of the electrodes and ensure that your skin is free of cream and fat. Do not clean your skin with alcohol.

⚠ **Warning:** Ensure that your iOS device is sufficiently charged in order to perform ECG readings.

⚠ **Warning:** During use of CardioSecur, do not connect any other device or computer via an adapter to your iOS device and/or to the headset jack of your iOS device.

⚠ **Warning:** During use of CardioSecur ensure that neither the patient nor any of the electrodes come into contact with other persons or conductive objects.

⚠ **Warning:** In case you cannot perform a firmware update, please contact Personal MedSystems' customer support.

⚠ **Warning:** Disinfect CardioSecur before each use. Please take account of warnings and information in section 12 Care and Maintenance,.

⚠ **Warning:** Perform regular data backups in order to prevent data loss (see section 5, Data Backup).

## 4 Getting Started

In order to run CardioSecur Pro on your iPad, download the free **CardioSecur Pro** app from the App-Store to your iPad:

## 5 Data Backup

You may perform a data backup in the following ways:

- Patient and ECG data may be stored and recovered together with the overall data from your iOS device, via iCloud or iTunes.
- ECG data may be exported in PDF format e. g. via e-mail to other databases (e. g. PMDS, HIS) or storage media (e. g. hard drive, memory card).

## 6 ISP

The ISP (integrated service platform) offers the following possibilities:

- Platform to connect physician and patient
- Save, store and transfer data between different devices
- Enable new features for the user


Different settings in the ISP are described in the following sections

- iPad:
- iPhone:

## 7 GDT-Interface

### 7.1 Overview

The GDT-interface can be used to transfer data between CardioSecur pro and an external data-base.

### 7.2 WebDAV-Server

The WebDav-Server needs directories with the following functions:

- Directory of jobs
- Directory of results

For this task there can be two separate directories or both functionalities can be placed in the same directory.

Setting up WebDav on Windows 7/8 Professional:

First, you must install the Internet Information Services (IIS) and WebDAV publishing features of Windows. In Windows 7, click Start, type "Windows features," and then click enter. In Windows 8, go to the Start screen, type "Windows features," select "Settings" on the right, and then select the shortcut that appears.

Make sure you select at least the following:

- Internet Information Services
- IIS Management Console
- WebDAV Publishing
- Windows Authentication

Next, to enable the WebDAV publishing you must open the IIS Manager. In Windows 7, click Start, type "IIS," and open the shortcut that appears. In Windows 8, go to the Start screen, type "IIS," and then select the shortcut that appears.

To configure the authentication, select the "Default Web Site" and double-click the "Authentication" icon. The list of authentication types displayed will depend upon which ones you installed earlier, but I'd suggest disabling "Anonymous" and enabling "Windows Authentication" if you will only be using IIS for WebDAV.

Next, select the "Default Web Site" site again and double-click the "WebDAV Authoring Rules" and then click "Enable" on the right. Now click "Add Authoring Rule" on the right to configure the permissions and make your desired selections.

To use https do a right click on the default website and add binding.

- Open IIS Manager and navigate to the level you want to manage. For information about opening IIS Manager, see [Open IIS Manager \(IIS 7\)](#). For information about navigating to locations in the UI, see [Navigation in IIS Manager \(IIS 7\)](#).
- In Features view, double-click Server Certificates.
- In the Actions pane, click Create Self-Signed Certificate.
- On the Create Self-Signed Certificate page, type a friendly name for the certificate in the Specify a friendly name for the certificate box, and then click OK.

Add WWW-Services to your firewall.

## 7.3 Data structures

### Job list

The name of the GDT-files for jobs have the following structure: "MOBDSxxx.yyy".

xxx and yyy are parameters that can be chosen independently to identify files from the patient database that are ignored by CardioSecur pro. Job files are deleted by CardioSecur pro after they have been read.

CardioSecur pro is able to use following files:

- Write ECG (GDT 6502)
- Master data (GDT 6501)

The content of the files of CardioSecur pro uses the parameters of the following table.

CardioSecur pro also uses the master data if it is included in the write ECG file.

Description	dataset identifier	Style/Content
dataset identifier	8000	alphanumeric
patient identifier	3000	alphanumeric
Patient first name	3102	alphanumeric
Patient last name	3101	alphanumeric
Patient date of birth	3103	style: DDMMYYYY DD = day MM = month YYYY = year z. B. 3.10.2016 -> 03102016

### Job Result

The name of the result GDT-file has the following structure: "DSMOBxxx"

xxx is a consecutive number to identify files that are provided by CardioSecur pro. Job requests are deleted by the receiver after they have been read.

Next to the GDT-file a PDF-report or a file with measurement data (SAV-file) in the SCP-format is created. The filename is included in the DSMOB.xxx file.

The result consists of two files:

- DSMOB.xxx with GDT-information
- binary files (PDF, SAV)

The content of the files is defined by the GDT-fields defined in the following table.

Description	Dataset identifier	Style/Content
Measurement identifier		„MOBEKG“
Device procedure identifier	8402	„EKG00“
Patient identifier	3000	= patient identifier of job list
Patient first name	3102	= patient first name from job
Patient last name	3101	= patient last name from job list
Date of the measurement	8432	format: DDMMYYYY DD = day MM = month YYYY = year z. B. 3.10.2016 -> 03102016
Time of the measurement	8439	format: HHMMSS HH = hour MM = minute SS = seconds e.g. 17:30:20s -> 170320
Archive identifier	6302	„1“
Format of the content	6303	<ul style="list-style-type: none"> <li>• „PDF“ for PDF-Report</li> <li>• „SAV“ for raw data in the SCP-format</li> </ul>
Description	6304	„mobile ECG“
File name Result file	6305	result file name

## 7.4 Process

- The host generates job files and stores these in the directory set up for the jobs
- CardioSecur pro reads the job files from the specified directory
- CardioSecur pro deletes the job file after the read process was successful
- These jobs are displayed in the patient database of the CardioSecur pro
- The physician writes ECGs for all patients that have a request in the patient database
- For completed jobs the GDT-file is stored in the directory specified for the results together with the PDF-file or the SAV-file.

## 8 Care and Maintenance

### Cleaning and Maintenance of iOS device

Take account of any information provided by Apple Inc. regarding cleaning and maintenance of your iOS device.

### Cleaning and disinfection of CardioSecur

Clean CardioSecur with a slightly dampened towel only. Cleaning with water and soap is usually sufficient. Any liquids on the casing need to be dried off with a towel immediately.

**⚠ Warning:** Never clean CardioSecur with benzene, benzine, paint thinner or concentrated alcohol.

To disinfect CardioSecur, wipe it with a disinfectant. Only use non-aggressive agents for disinfection. Recommended disinfectants (for surface disinfection as recommended by the respective manufacturers of these disinfectants in their user manuals) are Helipur® H plus N and Sagrotan® med concentrate.

**⚠ Warning:** Never dip CardioSecur or its parts into cleaning agents and do not sterilize it in a steam sterilizer (autoclave) or gas sterilizer (ethylene oxide, formaldehyde, ozone etc.).

### Maintenance of CardioSecur

**⚠ Warning:** Do not open CardioSecur. CardioSecur does not contain any parts that need to be maintained by the user.

## 9 Technical Data and Operating Conditions

Parameter	Description
Model	CardioSecur Pro
Lead Method	Bipolar, 3-lead
Reading quality	22-lead ECG with 4 electrodes
Pulse range	18 to 256 beats/minute
Precision	Bandwidth: 0,05 to 125 Hz, sampling rate: 250 Hz
Type of protection	IP22
Heart rate calculation	Evaluation of R wave across minimum of 3 beats, otherwise across 7 seconds
Filter	Selectable: muscle tension filter 35 Hz (-3 dB) and 50 Hz filter
Storage	External storage on iOS device. Each gigabyte capacity is 10,000 min ECG readings
Maximum duration of an ECG	1800 seconds = 30 minutes
Electrical safety	External power supply, application unit type BF
Power supply	External power supply via iOS device
Temperature and humidity operation	+5 °C to +40 °C, 10 % to 93 %, non-condensing
Temperature and humidity storage	-25 °C to +70 °C, 10 % to 93 %, non-condensing
Air pressure for operation	700 - 1060 hPa
Electrodes	Disposable electrodes

**⚠ Warning:** Do not store or operate CardioSecur under the following conditions:

- Direct sunlight
- Wet or damp places
- Dusty places
- Vicinity of fire or open flames
- Places subjected to strong physical shocks
- Places subjected to high electromagnetic radiation

**ℹ Information:** Safety-related controls are not necessary during the product's lifetime.

## 10 Electromagnetic Compliance (EMC)

**⚠ Warning:** Note important information on electromagnetic compliance (EMC). The growing number of electronic devices like PCs, displays, TV sets, washing machines and mobile phones may subject medical devices to electromagnetic interference. This may cause malfunctioning of the medical device and create potentially unsafe situations. Equally, medical devices should not interfere with other electronic devices. The norm IEC60601-1-2 was introduced to regulate the requirements for EMC and prevent unsafe product situations. This norm defines the degree of immunity against electromagnetic interference and the maximum electromagnetic emission levels for medical devices. CardioSecur complies with norm IEC60601-1-2:2001 concerning immunity as well as emissions. Nonetheless special precautions should be taken. Do not use devices that cause strong electromagnetic fields in the vicinity of CardioSecur. **A minimum distance of 2 meters is recommended. In case that the distance is less than 2 meters CardioSecur may not perform properly.** In accordance with the norm IEC60601-1-2:2001 you may contact our customer support (see section **Error! Reference source not found., Error! Reference source not found.**) for further documentation.

### DIN EN 550112

Industrial, scientific and medical high frequency devices (ISM devices)

Radio frequency interference - limit values and measuring procedure

Interference emission	Limit value class	Result
Interference voltage	--	n/a
RFI field strength (not accredited testing)	B	OK

**DIN EN 60601 1-1-2**

Medical electrical equipment

Part 1-2: General requirements for safety:

Collateral standard: electromagnetic compliance; requirements and testing

**DIN EN 60601 1-2-47**

Medical electrical equipment

Part 2-47: Specific requirements for safety including essential performance characteristics for outpatient electrocardiographic systems

Interference immunity against		Testing accuracy		Result
		normative	tested	
Static electrical discharge	air	8 kV	8 kV	OK
	contact	6 (8) kV	8 kV	OK
High frequency electromagnetic fields (80 - 2500 mHz)	housing	3 V/m	10 V/m	OK
Fast transient disturbances (chest)	mains	2 kV	2 kV	n/a
	I/O	1 kV	-	n/a
Surge voltages	unsym.	2 kV	-	n/a
	sym.	1 kV	1 kV	n/a
Conducted disturbances, induced by high frequency	mains	3 V	10 V	n/a
	I/O	3 V	-	n/a
Energy-related frequency - Magnetic field 50 Hz (Helmholtz coil)	housing	3 A/m	3 A/m	n/a
Voltage dips, short-term interruptions Test levels as % of UT duration in supply frequency cycles All testing at 100 V and 240 V each 50 Hz	mains UT 100 - 240V 50 / 60 Hz	<5 % 0,5 cycles	<5 % 0,5 cycles	n/a
		40 % 100 cycles	40 % 100 cycles	n/a
		70 % 500 cycles	70 % 500 cycles	n/a
		<5 % 5000 cycles	<5 % 5000 cycles	n/a

## 11 Disposal

⚠ **Warning:** Ensure proper disposal of CardioSecur (electronic waste).

Labelling on CardioSecur, its packaging, and this user manual indicate that CardioSecur must not be disposed together with ordinary household refuse after its lifetime. Dispose of CardioSecur separately, so that harm to the environment or human health is not caused by uncontrolled disposal.

## 12 Labelling

[www.cardiosecur.com](http://www.cardiosecur.com)  
Personal MedSystems GmbH  
Wilhelm-Leuschner-Str. 41  
60329 Frankfurt am Main  
Made in Germany 201x

## 13 Contact

**Personal MedSystems GmbH**  
Wilhelm-Leuschner-Str. 41  
60329 Frankfurt am Main  
Germany

Phone: +49-(0)69-9072013-0  
Fax: +49-(0)69-9072013-99  
E-mail: [info@cardiosecur.com](mailto:info@cardiosecur.com)  
Website: [www.mobile-ecg.com](http://www.mobile-ecg.com) / [www.cardiosecur.com](http://www.cardiosecur.com)

Refer to our customer support in case of queries regarding CardioSecur.

Phone: +49-(0)69-9072013-0  
E-mail: [info@cardiosecur.com](mailto:info@cardiosecur.com)

# 14 iPad

## 14.1 ECG Recording


The following flow diagram shows the steps on how to perform an ECG recording. Details to each step can be found in the chapters 14.4 Patient Administration, 14.5 Patient Details, 14.6 ECG Display and Functions and 14.8 Export.

1


Start the CardioSecur pro app.

2

Connect the CardioSecur pro ECG cable to your iPad.

3

Attach the disposable electrodes to the patient's body according to the diagram.

4


Create a new patient file by touching the icon  .

5


Enter the patient's personal data and start a reading by touching the icon  .

6

Start an ECG recording by touching the icon  .


7

Stop an ECG recording by touching the icon .

When a recording is stopped it is stored automatically in the archive.

8


Create a report by touching the icon .


You may also edit a report later by selecting the recording form the archive.

9


Complete your report entries by touching the icon .

10


Send an ECG report as PDF-file via e-mail by touching the icon .


11

Select the preferred display of the ECG recording and attach it to the e-mail by touching the icon .

12

Enter the recipient's e-mail address and start the transmission by touching the icon .


## 14.2 Emergency ECG Recording

In an emergency you have the possibility to perform an ECG recording fast and without need for entering the patient's personal data. This may save in some cases valuable time. To commence an emergency recording touch the icon  in the patient administration window:

CardioSecur *pro* thereupon creates by default a patient with the following properties:

- Last name: current date and time (format: YYYY-MM-DD - hh:mm:ss)
- First name: Emergency
- Gender: male

An ECG reading will commence instantly without need for entering further input. In order to perform an ECG recording in emergency cases, follow the instructions given in chapter 14.1 ECG Recording starting at step 6.

## 14.3 Menu


### Physician Data

Under this menu item you may enter all physician data. When creating an ECG report your data will be added automatically to the report.


Physician data comprises the following information:

- Title
- TitleFirst name
- Last name
- Phone number
- Organisation
- Department
- Street
- Street number
- ZIP code
- Town
- Country

Physician Data	
Title	Dr.
First name	Tiny
Last name	Dooittle
Phone number	+00-123456789
Organisation	Dooittle Clinic
Department	Cardiology
Street	Dooittle Street
Street number	1
ZIP code	11111
Town	Dooittle Town
Country	Fairyland

<p>ISP Settings</p>	<p>Connecting to your personal user account on the ISP (Integratet Service Platform) allows you to use all your CardioSecur features. The ISP also is the center point of your connection to patients that use the CardioSecur active.</p>	
<p>Automatic ECG Interpretation</p>	<p>The automatic interpretation is a complete diagnostic tool that enables you to perform a full automated diagnostic analysis on a 10 second measurement. The interpretation also holds measurement values for characteristic ECG parameters as QRS duration and QT-time.</p>	
<p>GDT worklist</p>	<p>The GDT interface holds the possibility to integrate CardioSecur in your IT-System. Patients from your database can be synchronized into your CardioSecur app. ECGs taken for these patients can then be synchronized to your patient files of your IT system as PDF or as raw data.</p>	


<p><b>Intended Use</b></p>	<p>You will find the Intended Use in the app as well as in this user manual: see Chapter 2 Intended Use</p>	
<p><b>Safety Information</b></p>	<p>You will find the safety information in the app as well as in this user manual: ➤ See chapter 3 Safety Information</p>	
<p><b>About CardioSecure pro</b></p>	<p>Information about CardioSecure pro comprises the following Data:</p> <ul style="list-style-type: none"> <li>• Serial number</li> <li>• Version of the application</li> <li>• Version of the ECG cable</li> <li>• Version of the firmware</li> <li>• Contact information of Personal MedSystems</li> <li>• Labelling</li> </ul>	

## 14.4 Patient Administration

Menu	➤ See chapter 14.3 Menu
	Performing an emergency recording ➤ See chapter 14.2 Emergency ECG Recording
Touch a patient entry in the list	Access a patient's details ➤ See chapter 15.5 Patient Details
Add Patient	Create a new patient ➤ See chapter 15.5 Patient Details
With one finger swipe from right to left across the patient entry (delete gesture)	Delete a patient entry by touching the icon  .  
Symboles	
	The ISP-Symbol symbolize patients that are synchronized with the ISP (See chapter 14.9 Patient Details)


## 14.5 Patient Details

<div style="display: flex; justify-content: space-between; align-items: center;"> <span>No SIM</span> <span>99%</span> </div> <div style="display: flex; justify-content: space-between; align-items: center; background-color: #e0f2f1; padding: 5px;"> <span>Patients</span> <span>John Doe</span> </div>							
<div style="display: flex; justify-content: space-around; margin-bottom: 10px;"> <div style="border: 1px solid #ccc; padding: 10px; background-color: #e0f2f1; display: flex; align-items: center; gap: 10px;">  <span>New reading</span> </div> <div style="border: 1px solid #ccc; padding: 10px; background-color: #e0f2f1; display: flex; align-items: center; gap: 10px;">  <span>Archive</span> </div> </div> <div style="background-color: #e0f2f1; padding: 10px; margin-bottom: 10px;"> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 20%; padding: 5px;">Patient ID</td> <td style="padding: 5px;">0123456789</td> </tr> <tr> <td style="padding: 5px;">Last name</td> <td style="padding: 5px;">Doe</td> </tr> <tr> <td style="padding: 5px;">First name</td> <td style="padding: 5px;">John</td> </tr> </table> </div>		Patient ID	0123456789	Last name	Doe	First name	John
Patient ID	0123456789						
Last name	Doe						
First name	John						
<b>New reading</b>	Commence a new ECG reading for a patient ➤ See chapter 14.6 ECG Display and Functions						
<b>Archive</b>	List of all ECG recordings of a patient						
<b>Data fields</b>	<p>You may enter personal data of the patient. The patient's data will be added to each performed ECG recording.</p> <p>Data concerning the patient comprises the following information:</p> <ul style="list-style-type: none"> <li>• Patient ID</li> <li>• Last name</li> <li>• First name</li> <li>• Date of birth</li> <li>• Height</li> <li>• Weight</li> <li>• Gender</li> <li>• E-mail address</li> <li>• Mobile phone number</li> <li>• Phone number</li> <li>• Fax number</li> <li>• Street</li> <li>• Street number</li> <li>• ZIP code</li> <li>• Town</li> <li>• Country</li> </ul>						

## 14.6 ECG Display and Functions

ECG display may be performed in the following modes:

- **Viewing mode:** In viewing mode ECG readings can be viewed. You have the option to halt the running ECG reading and resume it again. ECG readings are merely displayed not recorded in this mode.
- **Recording mode:** In recording mode ECG readings can be recorded and stored in the local memory of your iPad.


### Viewing mode


	<p>Halt a running ECG reading. This icon appears, while the ECG reading is running (not in recording mode).</p>
	<p>Resume a running ECG reading. This icon appears, when a running ECG image was halted.</p>
	<p>Start an ECG recording. This icon appears after activating a new reading.</p>

## Recording mode

Stop an ECG recording. This icon appears, when an ECG recording was commenced. By stopping an ECG recording it will be stored in the archive.

## Report Function

Either after performing an ECG recording or when selecting a recording from the archive you may create a report attaching to the recording.

Reports allow you to edit the following options:

- Rhythm
- QRS complex
- Repolarisation
- Summary
- Comment


## General Functions

Display options:

- Number of columns
- Lead block selection
- Individual leads
- Cabrera display

<b>Time axis</b>	Selection of time axis scaling: <ul style="list-style-type: none"> <li>• 5 mm/s</li> <li>• 10 mm/s</li> <li>• 25 mm/s</li> <li>• 50 mm/s</li> </ul>
<b>Voltage axis</b>	Selection of voltage axis scaling (amplitude): <ul style="list-style-type: none"> <li>• 0.5 mV/cm</li> <li>• 1 mV/cm</li> <li>• 2 mV/cm</li> <li>• 5mV/mm</li> </ul>
<b>Filter</b>	Selection of ECG filters: <ul style="list-style-type: none"> <li>• No filter</li> <li>• 50 Hz filter</li> <li>• Muscle tension filter</li> </ul>
Swipe with one finger (scrolling gesture)	Scrolling on the time axis. The header displays the current view position in the ECG recording (current/total duration in s). 
Narrow and widen space between two fingers (zoom gesture)	Zooming in and out the current view.



## 14.7 Archive

<p>Touch an entry</p>	<p>Activate the display of an archived ECG recording.          ➤ See chapter 14.6 ECG Display and Functions.</p>
<p>Swipe with one finger from right to left across an entry (delete gesture)</p>	<p>Delete an archived ECG recording with the icon .</p>

## 14.8 Export

Every ECG may be exported via e-mail after recording or from the archive. In order to do so an e-mail account must be installed on your iPad. The ECG report attaches to the e-mail as PDF-file. The alignment of an ECG report is by default and unalterable in landscape view. Follow the subsequent steps to export an ECG recording:

Follow the subsequent steps to export an ECG recording:




<p><b>Export</b></p>	<p>Every ECG recording that is new or selected from the archive may be exported. To do so touch the <b>Export</b> icon in the header.</p>  <p>The screenshot shows the top header of the application. It includes a status bar at the top with 'No SIM' and '99%' battery. Below that is a navigation bar with 'Archive' on the left, the date and time '2012-12-23 - 14:15:24 - 14s/16s' in the center, and an 'Export' button on the right. Below the navigation bar, the patient's name 'John Doe' and vital signs '01.01.1951   190cm   95kg' are displayed.</p>
<p><b>E-mail</b></p>	<p>Before exporting you may select how the ECG recording is displayed in the PDF-file.</p> <p>Displaying options:</p> <ul style="list-style-type: none"> <li>• Number of columns</li> <li>• Lead block selection</li> <li>• Individual leads</li> <li>• Cabrera display</li> <li>• Rhythm strip</li> <li>• Entire ECG or current view</li> </ul> <p>Complete your selection by touching the <b>E-Mail</b> icon and the e-mail send dialog will appear.</p>  <p>The screenshot shows the 'Leads' selection screen. It has two columns: 'Leads' and 'Columns'. Under 'Leads', there are checkboxes for various lead blocks: V1-V3, V4-V6, aVR, aVL, aVF, I, II, III, R, RL, L, V5-V6, V7-V9, Y1-Y6, V5I-V6I, and aVL, aVF, aVF, I, II, III. Under 'Columns', there are checkboxes for 1, 2, 3, 4, 5, 6, 7, 8, 9, and 10 columns. There is also a 'Cabrera' option at the bottom.</p>
<p><b>Send</b></p>	<p>The ECG report will be attached to the e-mail as PDF-file.</p> <p>The e-mail subject will contain the patient's name. You may change the subject and content of the e-mail. In order to send the e-mail you need to enter the e-mail address of the intended recipient.</p> <p>Send the e-mail by touching the icon <b>Send</b> the e-mail will be sent.</p>  <p>The screenshot shows an e-mail composition screen. The 'To:' field is filled with 'andreas@supercard.com'. The 'Subject:' field is filled with 'John Doe'. There is a 'Send' button at the bottom right.</p>

# 15 iPhone

## 15.1 ECG Recording


The following flow diagram shows the steps on how to perform an ECG recording. Details to each step can be found in the chapters 14.4 Patient Administration, 14.5 Patient Details, 14.6 ECG Display and Functions and 14.8 Export.

<p>1</p>  <p>Start the CardioSecur pro app.</p>	<p>2</p>  <p>Connect the CardioSecur pro ECG cable to your iPad.</p>	<p>3</p>  <p>Attach the disposable electrodes to the patient's body according to the diagram.</p>
<p>4</p>  <p>Create a new patient file by touching the icon .</p>	<p>5</p>  <p>Enter the patient's personal data and start a reading by touching the icon .</p>	<p>6</p>  <p>Start an ECG recording by touching the icon .</p>


Stop an ECG recording by touching the icon .

When a recording is stopped it is stored automatically in the archive.


Create a report by touching the icon .

You may also edit a report later by selecting the recording from the archive.


Complete your report entries by touching the icon .


Reading

Send an ECG report as PDF-file via e-mail by touching the icon .



Select the preferred display of the ECG recording and attach it to the e-mail by touching the icon .

Enter the recipient's e-mail address and start the transmission by touching the icon **Send**.

## 15.2 Quick reading

In an emergency you have the possibility to perform an ECG recording fast and without need for entering the patient's personal data. This may save in some cases valuable time. To commence an emergency recording touch the icon  in the patient administration window:



CardioSecur *pro* thereupon creates by default a patient with the following properties:



- Last name: current date and time (format: YYYY-MM-DD - hh:mm:ss)
- First name: Emergency
- Gender: male



An ECG reading will commence instantly without need for entering further input. In order to perform an ECG recording in emergency cases, follow the instructions given in chapter 14.1 ECG Recording starting at step 6.

## 15.3 Menu



		
<p><b>Physician Data</b></p>	<p>Under this menu item you may enter all physician data. When creating an ECG report your data will be added automatically to the report.</p> <p>Physician data comprises the following information:</p> <ul style="list-style-type: none"><li>• Title</li><li>• First name</li><li>• Last name</li><li>• Phone number</li><li>• Mobile number</li><li>• Fax number</li><li>• E-mail</li><li>• Organisation</li><li>• Department</li><li>• Street</li><li>• Street number</li><li>• ZIP code</li><li>• Town</li><li>• Country</li></ul>	

<p><b>ISP Settings</b></p>	<p>Entering the personal data and synchronizing will allow you to use ISP. Using the ISP features can be added or removed, and ECG data generated by patients can be viewed.</p> 
<p><b>Automatic ECG interpretation</b></p>	<p>This allows you to activate/deactivate the automatic interpretation, which is a diagnostic ECG program of the analysis of different ECG data.</p> 

<p><b>GDT Worklist</b></p>	<p>The GDT interface holds the possibility to integrate CardioSecur in your IT-System. Patients from your database can be synchronized into your CardioSecur app. ECGs taken for these patients can then be synchronized to your patient files of your IT system as PDF or as raw data.</p>	
<p><b>Intended Use</b></p>	<p>You will find the Intended Use in the app as well as in this user manual: see Chapter 2</p>	 <p>CardioSecur is intended to record, evaluate, document and diagnose a resting ECG of adults.</p> <p>CardioSecur is a medical, electrical system comprising an ECG cable, an iPad and an app for an iPad. CardioSecur records a 12-lead ECG with four electrodes.</p> <p>Only medically trained personnel may use CardioSecur.</p> <p>In order to record an ECG with CardioSecur, the disposable electrodes need to be positioned on the patient's thorax according to the following diagram:</p> <ul style="list-style-type: none"> <li>White electrode On the upper end of the sternum</li> <li>Yellow electrode On the lower end of sternum</li> <li>Red electrode In the middle of the right shoulder at the level of the lower end of the sternum</li> </ul>

<p><b>Safety Information</b></p>	<p>You will find the safety information in the app as well as in this user manual:</p> <p>➤ See chapter 3 Safety Information</p>	
<p><b>About CardioSecur pro</b></p>	<p>Information about CardioSecur <i>pro</i> comprises the following Data:</p> <ul style="list-style-type: none"> <li>• Serial number</li> <li>• Version of the application</li> <li>• Version of the ECG cable</li> <li>• Version of the firmware</li> <li>• Contact information of Personal MedSystems</li> <li>• Labelling</li> </ul>	

## 15.4 Patient Administration

	
<p><b>Menu</b></p>	<p>➤ See chapter 14.3 Menu</p>
	<p>Performing an emergency recording</p> <p>➤ See chapter 14.2 Emergency ECG Recording</p>
<p>Touch a patient entry in the list</p>	<p>Access a patient's details</p> <p>➤ See chapter 14.5 Patient Details</p>
	<p>Create a new patient</p> <p>➤ See chapter 14.5 Patient Details</p>
<p>With one finger swipe from right to left across the patient entry (delete gesture)</p>	<p>Delete a patient entry by touching the icon </p> 


## 15.5 Patient Details

	
	<p>Commence a new ECG reading for a patient</p> <p>➤ See chapter 14.6 ECG Display and Functions</p>
	<p>List of all ECG recordings of a patient</p>
<p>Data fields</p>	<p>You may enter personal data of the patient. The patient's data will be added to each performed ECG recording.</p> <p>Data concerning the patient comprises the following information:</p> <ul style="list-style-type: none"><li>• Patient ID</li><li>• Last name</li><li>• First name</li><li>• Date of birth</li><li>• Height</li><li>• Weight</li><li>• Gender</li><li>• E-mail address</li><li>• Mobile phone number</li><li>• Phone number</li><li>• Fax number</li><li>• Street</li><li>• Street number</li><li>• ZIP code</li><li>• Town</li><li>• Country</li></ul>


## 15.6 ECG Display and Functions

ECG display may be performed in the following modes:

- **Viewing mode:** In viewing mode ECG readings can be viewed. You have the option to halt the running ECG reading and resume it again. ECG readings are merely displayed not recorded in this mode.
- **Recording mode:** In recording mode ECG readings can be recorded and stored in the local memory of your iPad.

### Viewing mode

Halt a running ECG reading. This icon appears, while the ECG reading is running (not in recording mode).

Resume a running ECG reading. This icon appears, when a running ECG image was halted.

Start an ECG recording. This icon appears after activating a new reading.

## Recording mode

Axis/Filter

Leads

Stop Recording

Stop an ECG recording. This icon appears, when an ECG recording was commenced. By stopping an ECG recording it will be stored in the archive.

## Report Function

Axis/Filter

Leads


Report


Either after performing an ECG recording or when selecting a recording from the archive you may create a report attaching to the recording.

Reports allow you to edit the following options:

- Rhythm
- QRS complex
- Repolarisation
- Summary
- Remarks


# Report Function



## Display options:

- Number of columns
- Lead block selection
- Individual leads
- Cabrera display

Axis/Filter		
   		
<p><b>Time axis</b></p>	<p>Selection of time axis scaling:</p> <ul style="list-style-type: none"> <li>• 5 mm/s</li> <li>• 10 mm/s</li> <li>• 25 mm/s</li> <li>• 50 mm/s</li> </ul>	
<p><b>Voltage axis</b></p>	<p>Selection of voltage axis scaling (amplitude):</p> <ul style="list-style-type: none"> <li>• 0.5 mV/cm</li> <li>• 1 mV/cm</li> <li>• 2 mV/cm</li> <li>• 5mV/cm</li> </ul>	
<p><b>Filter</b></p>	<p>Selection of ECG filters:</p> <ul style="list-style-type: none"> <li>• No filter</li> <li>• 50 Hz filter</li> <li>• Muscle tension filter</li> </ul>	
General Functions		
<p>Swipe with one finger (scrolling gesture)</p>	<p>Scrolling on the time axis. The header displays the current view position in the ECG recording (current/total duration in s).</p> 	
<p>Narrow and widen space between two fingers (zoom gesture)</p>	<p>Zooming in and out the current view.</p>	



## 15.7 Archive

	
Touch an entry	Activate the display of an archived ECG recording. ➤ See chapter 15.6 ECG Display and Functions.
Swipe with one finger from right to left across an entry (delete gesture)	Delete an archived ECG recording with the icon  . 

## 15.8 Export

Every ECG may be exported via e-mail after recording or from the archive. In order to do so an e-mail account must be installed on your iPhone. The ECG report attaches to the e-mail as PDF-file. The alignment of an ECG report is by default and unalterable in landscape view.

Follow the subsequent steps to export an ECG recording:


<b>Export</b>	<p>Every ECG recording that is new or selected from the archive may be exported. To do so touch the  icon in the footer.</p>  <p>The footer contains four icons: a square with an arrow pointing right (Axis/Filter), two interlocking gears (Leads), a square with an arrow pointing right (Export), and a clipboard with a document (Report).</p>
<b>E-mail</b>	<p>Before exporting you may select how the ECG recording is displayed in the PDF-file. 15.8.1.1.1.1</p> <p>Displaying options:</p> <ul style="list-style-type: none"><li>• Number of columns</li><li>• Lead block selection</li><li>• Individual leads</li><li>• Cabrera display</li><li>• Rhythm strip</li><li>• Entire ECG or current view</li></ul> <p>Complete your selection by touching the  icon and the e-mail send dialog will appear.</p>  <p>The screenshot shows the 'Export' menu with the following options: Leads (V1 - V6 aVL aVR aVF I III), Columns (2), Cabrera (No), Rhythm strip (No rhythm strip), and Area (Entire ECG).</p>

**Send**

The ECG report will be attached to the e-mail as PDF-file.

The e-mail subject will contain the patient's name. You may change the subject and content of the e-mail. In order to send the e-mail you need to enter the e-mail address of the intended recipient.

Send the e-mail by touching the icon  the e-mail will be sent.

[www. cardiosecur.com](http://www.cardiosecur.com)  
Personal MedSystems GmbH  
Wilhelm Leuschner Straße 41  
60329 Frankfurt am Main  
Deutschland/Germany