CardioSecur *ACTIVE*

Knows your Heart

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User Manual (EN)

Version V5.0

17.12.2015

1 Index

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Telefon: +49-(0)69-907477-81
 E-Mail: info@cardiosecur.com

1 Index

1	Index	1
2	Symbols	2
3	Intended Use	3
4	Safety Information	5
5	Connecting the CardioSecur <i>ACTIVE</i> ECG-cable	8
6	Disposable Electrodes - P-00-S/xx	9
7	Serial Number and Contact Address	11
8	Care and Maintenance	12
9	Technical Data and Operating Conditions	13
10	Electromagnetic Compliance - EMC	14
11	Disposal	16
12	Labelling	16
13	Contact	16

2 Symbols

Packaging or instructions included in the scope of delivery may display the following symbols:

Symbol	Meaning
	Use disposable electrodes before the displayed date of expiry.
LOT	Lot denomination
	This product complies with the requirements of the directive for medical devices 93/42/EWG.
REF	Item number
	Attention Attention concerns information of high relevance.
	Warning Warning concerns safety-relevant information.
	Information Information concerns details that are worth knowing.
	Manufacturer's data
	Store at indicated temperature with unbroken seal.
	Only for single use
	Latex-free
	CardioSecur <i>ACTIVE</i> and its accessories - like any other electronic device - contain metal and synthetic parts, which need to be disposed separately after operating life in order to prevent negative environmental impact. Send CardioSecur <i>ACTIVE</i> and its accessories with sufficient postage and if possible in its original packaging to Personal MedSystems GmbH for free of charge disposal.
	Application unit of type BF. An application unit of type BF grants protection against electric shock, in particular with regard to permissible discharge current, according to the requirements of norm EN 60601-1.

3 Intended Use


CardioSecur is intended to record, evaluate, document and diagnose a resting ECG of adults.

CardioSecur is a medical electrical system, comprising an ECG cable, an App for Smartphone, iPod touch or iPad. CardioSecur can record, evaluate and store a complete 12-lead ECG with four electrodes. CardioSecur is able to detect arrhythmia and ischemia of the heart by comparing a reference reading with a control reading. CardioSecur serves the purpose of self-monitoring heart activity.

The user is required to record a personal reference reading in a state free of health problems and discomfort. The reference reading is stored in CardioSecur. The user may perform the reading himself or under medical supervision. A control reading may be performed either immediately in case of discomforting symptoms or in regular intervals (e.g. every other day) independently of discomforting symptoms.

The user is not required to have previous medical knowledge or special training for using CardioSecur. CardioSecur is recommended to individuals with a known or suspected heart condition or bearing risk factors for cardiac disease (e.g. hypertension, diabetes mellitus or increased cholesterol levels). Patients and individuals taking precautionary measures may use CardioSecur.

In order to record an ECG (reference and control reading), electrodes need to be positioned on the user's thorax according to the following diagram:

Prior to a reading the user is required to sit calmly and relax for a minimum of 3 minutes. During a reading (duration ca. 1 minute) the user is continuously required to sit calmly and relax without talking or coughing.

CardioSecur *ACTIVE* performs a comparative evaluation between the reference reading and a control reading and thereupon provides the user a recommendation to act:

- **White result:** resembles “no changes” compared to the personal reference ECG. There is generally no need for action.
- **Yellow result:** resembles “minor irregularities” compared to personal reference ECG. CardioSecur *ACTIVE* recommends planning for a doctor visit for clarification.
- **Red result:** resembles “severe irregularities” compared to reference ECG. CardioSecur *ACTIVE* recommends seeing a doctor instantly.

In accordance with applicable guidelines of the European Society of Cardiology, CardioSecur *ACTIVE* recommends visiting a doctor instantly whenever symptoms or discomfort prevail for more than 20 minutes. This applies despite the result given by CardioSecur *ACTIVE* shows “no changes” or “minor irregularities” compared to the reference ECG.

The user may store, print and administer all ECG readings (reference and control readings) taken with CardioSecur *ACTIVE* by means of a personal user account on the CardioSecur database, an optional service by Personal MedSystems. In addition users may self-select a physician to view personal ECG data for informational purposes.

This manual describes usage of CardioSecur *ACTIVE* on iPhone 4S, 5/5S/5C, 6/6plus (iPhone is a registered trademark of Apple Inc.) as well as for chosen Android smartphones. The current list can be found on www.cardiosecur.com. For reasons of simplicity we refer herein exclusively to smartphone.

4 Safety Information

Take account of the following warnings and information. In case of your non-compliance you run the risk of causing malfunction, damage and injury.

⚠ Warning: Only use CardioSecur *ACTIVE* in accordance with its intended use.

⚠ Warning: CardioSecur *ACTIVE* contains personalized data. It must only be used by the individual whose personal reference ECG is stored in it.

⚠ Warning: The results provided by CardioSecur *ACTIVE* must not be used as basis for the beginning or changing of therapeutic measures without an independent medical examination.

⚠ Warning: ECG readings by CardioSecur *ACTIVE* are not able to detect all health changes of cardiac disease! If you notice changes to your state of health, contact your physician instantly. In case of prevailing symptoms for more than 20 minutes, contact your doctor, even if the results provided by CardioSecur *ACTIVE* show no or minor irregularities compared to your reference ECG.

⚠ Warning: Do not attempt to make a diagnosis on the basis of the results and analysis provided by CardioSecur *ACTIVE*. Always check with your physician. Self-diagnosis may cause a deterioration of your state of health or medical condition.

⚠ Warning: Provided results by CardioSecur *ACTIVE* reflect the very moment of a single recording. Your state of health may change rapidly. In case you notice a change in your state of health contact a physician.

⚠ Warning: CardioSecur *ACTIVE* serves the purpose of detecting newly occurred ischemia and arrhythmia of the heart. CardioSecur *ACTIVE* is not designed for continuous monitoring of your ECG, but evaluates merely a short time window of your heart activity.

⚠ Warning: Note that a physician may only receive your ECG recordings and your attached comments if you are actively linked to the physician through the CardioSecur database. Transmission of ECG data to the linked physician is not guaranteed. In addition, when your ECG data is received by your physician, this does not guarantee that he or she will perform a diagnosis or look at the data. It is recommended that you contact him or her directly. This is especially important in the case of an emergency.

⚠ Warning: Before a reading with CardioSecur *ACTIVE* sit calmly for a minimum of 3 minutes. Ensure that your posture is relaxed. Being seated with a straight back is ideal. Do not cough or speak during a reading. Body movements may cause interference influencing the evaluation.

⚠ Warning: Do not perform readings when subjected to concussions, e.g. in a moving vehicle, train or plane.

⚠ Warning: Do not perform ECG readings if the attachment points of the electrodes on your skin are wet, e.g. due to sweat or after bathing.

⚠ Warning: Proper functioning of the disposable electrodes is influenced by intense body hair. Remove regularly hair from your skin at the attachment points of the electrodes.

⚠ **Warning:** Do not record a reference ECG with incorrectly positioned or transposed electrodes, or else a correct evaluation of a subsequent control reading is not possible.

⚠ **Warning:** You must renew your reference ECG after every acute event or intervention (e.g. heart attack, cardioversion in case of arrhythmia of the heart, heart bypass intervention, heart pace maker implantation, cardiac catheter intervention, heart valve intervention) or when recommended by your physician.

⚠ **Warning:** Only use the disposable electrodes type as recommended by Personal MedSystems (P-00S/xx).

⚠ **Warning:** Only use disposable electrodes that are unused, have not expired and are at room temperature, in order to ensure sufficient signal quality and prevent distorted results.

⚠ **Warning:** Remove all disposable electrodes after a reading from your skin, as prolonged adhesion may cause skin irritation.

⚠ **Warning:** Do not use CardioSecur *ACTIVE* in case of sensitive skin or allergies. This could cause rash, irritation of the skin or eczemas.

⚠ **Warning:** Only use the CardioSecur *ACTIVE* ECG cable as provided by Personal MedSystems.

⚠ **Warning:** Do not use CardioSecur *ACTIVE* in combination with external defibrillators. CardioSecur *ACTIVE* is not defibrillator proof.

⚠ **Warning:** During use of CardioSecur *ACTIVE*, do not connect any other device or computer via an adapter to your smartphone and/or to the headset jack of your smartphone.

⚠ **Warning:** Do not perform readings in places exposing CardioSecur *ACTIVE* to high electromagnetic radiation (e.g. TV-sets, monitors, dish washers).

⚠ **Warning:** Do not expose CardioSecur *ACTIVE* to any electrostatics. Discharge any electrostatics out of your body before touching CardioSecur *ACTIVE*.

⚠ **Warning:** Do not expose CardioSecur *ACTIVE* to any physical shocks or vibrations. Do not drop it on the floor or step on it.

⚠ **Warning:** Do not place anything onto CardioSecur *ACTIVE* or it might get damaged.

⚠ **Warning:** Do not open, take apart or repair CardioSecur *ACTIVE*. Do not transform CardioSecur *ACTIVE*. CardioSecur *ACTIVE* does not contain any parts that need to be maintained by the user.

⚠ **Warning:** Store and use CardioSecur *ACTIVE* only under the conditions described under section Technical Data.

⚠ **Warning:** Do not sterilize CardioSecur *ACTIVE* in a steam sterilizer (autoclave) or in a gas sterilizer (ethylene oxide, formaldehyde, ozone etc.)

⚠ **Warning:** Do not use CardioSecur *ACTIVE* in the vicinity of combustible anaesthetics, drugs or pressurized oxygen.

⚠ **Warning:** Only use CardioSecur *ACTIVE* if the ECG cable has no apparent damage.

⚠ **Warning:** Note the information in section 8, Care and Maintenance.




⚠ **Warning:** CardioSecur *ACTIVE* utilizes an ECG method computing a 12-lead ECG from 3 leads read with 4 electrodes. Deviations might arise in comparison to a 12-lead ECG read with 10 electrodes.

⚠ **Attention:** Use CardioSecur *ACTIVE* only with correctly set time and date.

⚠ **Attention:** Contact our customer support if you cannot manage to perform a firmware update.

5 Connecting the CardioSecur ACTIVE ECG-cable

Connect the CardioSecur ACTIVE ECG-cable to your smartphone.

			
<p>A CardioSecur ACTIVE with 30-Pin connector can be plugged directly into an iPhone with 30-Pin receptacle (e.g. iPhone 4).</p>	<p>A CardioSecur ACTIVE with 30-Pin connector can be plugged into an iPhone with Lightning receptacle (e.g. iPhone 5) only by means of the Apple Lightning adapter.</p> <p>This adapter is not part of the delivery contents!</p>	<p>A CardioSecur ACTIVE with Lightning connector can be plugged directly into an iPhone with Lightning receptacle (e.g. iPhone 5).</p>	<p>A CardioSecur ACTIVE with Micro-USB connector can be plugged directly into an Android smartphone (e.g. Samsung Galaxy S5).</p>

6 Disposable Electrodes - P-00-S/xx

Disposable electrodes constitute the link between your body and CardioSecur *ACTIVE*. They are of critical importance for the technology's proper function. Personal MedSystems cooperates with highly specialized manufacturers for disposable electrodes, in order to ensure optimal signal quality.

Intended Use

P-00-S/xx are disposable electrodes designed for reading electrocardiograms (ECG).

⚠ **Warning:** Use disposable electrodes only if the expiry date has not passed.


⚠ **Warning:** Use disposable electrodes only if they have reached room temperature.

⚠ **Warning:** Store disposable electrodes only with unbroken seal and at temperatures between 10°C and 25°C.

Operating Instructions


Use disposable electrodes P-00-S/xx according to the subsequent steps:

1

First connect the four ECG electrode connectors with each an unused disposable electrode. Align the protruding blue connector in direction of the ECG cable.

2

Remove the protective foil from the first disposable electrode and attach it without touching the adhesive side and without putting pressure on the central gel pad. Attach one by one the second, third and fourth disposable electrode in the same way and according to the diagram on your thorax (see section 3, Intended Use).

3

After a reading remove the disposable electrodes carefully from your skin. Careless removal may cause damage to your skin. Should you encounter difficulties in removing the disposable electrodes, dissolve the adhesive agent of the electrodes with water. Remaining gel may be removed from your skin with paper tissues or water.

Note the instructions provided on the backside of the disposable electrodes' pouch.

Manufacturer Information

Ambu A/S
Baltorpbakken 13
DK-2750 Ballerup
Denmark



Made in Malaysia

P-00-S/xx

7 Serial Number and Contact Address

NE

3

The following information is displayed on the Info page:

- Software version of your iPhone/Android app
- Hard- and software version of your CardioSecur *ACTIVE* ECG-cable
- Serial number of your CardioSecur *ACTIVE* ECG-cable
- Contact details of Personal MedSystems GmbH
- Marking of the product

8 Care and Maintenance

Cleaning and Maintenance smartphone

Attend the instructions by Apple Inc. in order to clean and maintain your smartphone.

Cleaning CardioSecur *ACTIVE*

Clean CardioSecur *ACTIVE* ECG-cable with a slightly dampened towel. In most cases cleaning with water and soap will be sufficient. Ensure that no liquids whatsoever penetrate the casing. Liquids on the casing must instantly be wiped off with a dry towel.

⚠ **Warning:** Clean CardioSecur *ACTIVE* never with benzene, benzine, paint thinner or concentrated alcohol.

Maintenance CardioSecur *ACTIVE*

⚠ **Warning:** Do not open CardioSecur *ACTIVE*. CardioSecur *ACTIVE* does not contain any parts that need to be maintained by the user.

Disinfection CardioSecur *ACTIVE*

For disinfecting the ECG-cable wipe it off with a disinfectant. For disinfecting CardioSecur *ACTIVE* use a non-aggressive solution. Recommended disinfectants (for superficial disinfection as recommended in the respective manufacturer's user manual) are:

- Helipur® H plus N
- Sagrotan® med Konzentrat

⚠ **Warning:** Never dip CardioSecur *ACTIVE* or parts of it into a cleaning liquid und do not sterilize CardioSecur *ACTIVE* in a steam sterilizer (autoclave) or in a gas sterilizer (ethylene oxide, formaldehyde, ozone etc.).

9 Technical Data and Operating Conditions

Parameter	Description
Model	CardioSecur <i>ACTIVE</i>
Lead method	Bipolar, 3 lead
Reading quality	12-lead ECG with 4 electrodes
Pulse range	18 to 256 beats/minute
Precision	Band width: 0,05 to 40 Hz, sampling rate: 250 Hz
Type of protection	IP20
Storage	10.000 readings
Electrical safety	Application unit type BF
Power supply	through Smartphone/iPod
Temperature and humidity for operation	+5 to +45 °C, 10 to 95 %, non-condensing
Temperature and humidity for storage	-20 to +45 °C, 10 to 95 %, non-condensing
Air pressure for operation	700 - 1060 hPa
Electrodes	Disposable electrodes, type P-00-S/xx

Warning: Do not store CardioSecur *ACTIVE* ECG-cable under the subsequent conditions:

- Direct sunlight
- Wet or damp places
- Dusty places
- Vicinity of fire or open flames
- Places subject to physical shocks
- Places subject to high electromagnetic radiation

Information:

- Keep the original packaging of CardioSecur *ACTIVE* and its accessories in order to identify the product in case of queries.
- We reserve the right to change technical data without prior notice.
- During product lifetime no safety-related controls are necessary.

10 Electromagnetic Compliance - EMC

⚠ **Warning:** Note important information on electromagnetic compliance (EMC). The growing number of electronic devices like PCs, displays, TV-sets and washing machines may subject medical devices to electromagnetic interference. This may cause malfunctioning of the medical device and create potentially unsafe situations. Equally medical devices should not interfere with other electronic devices. The norm IEC60601-1-2 was introduced to regulate the requirements for EMC and prevent unsafe product situations. This norm defines the degrees of immunity against electromagnetic interference and the maximum electromagnetic emission levels for medical devices. CardioSecur *ACTIVE* complies with norm IEC60601-1-2:2001 concerning immunity as well as emissions. Nonetheless special precautions should be taken: Do not use devices that cause strong electric or electromagnetic fields in the vicinity of CardioSecur *ACTIVE*. A minimum distance of 2 meters is recommended. In case that the distance is less than 2 meters CardioSecur *ACTIVE* may not perform properly. In accordance with the norm IEC60601-1-2:2001 you may request our customer support (see section 13, Contact) for further documentation.

DIN EN 550112

Industrial, scientific and medical high frequency devices (ISM devices)
Radio frequency interference - limit values and measuring procedure

Interference emission	Limit value class	Result
Interference voltage	--	n/a
RFI field strength (not accredited testing)	B	OK

DIN EN 60601 1-1-2

Medical electrical equipment

Part 1-2: General requirements for safety:

Collateral standard: electromagnetic compliance; requirements and testing

DIN EN 60601 1-2-47

Medical electrical equipment

Part 2-47: Specific requirements for safety including essential performance characteristics for outpatient electrocardiographic systems

Interference immunity against		Testing accuracy		Result
		normative	tested	
Static electrical discharge	air	8 kV	8kV	OK
	contact	6 (8) kV	8 kV	OK
High frequency electromagnetic fields (80 - 2500 mHz)	housing	3 V/m	10 V/m	OK
Fast transient disturbances (chest)	mains	2 kV	2 kV	OK
	I/O	1 kV	-	OK
Surge voltages	unsym.	2 kV	-	OK
	sym.	1 kV	1kV	OK
Conducted disturbances, induced by high frequency	mains	3 V	10 V	OK
	I/O	3 V	-	OK
Energy-related frequency - Magnetic field 50 Hz (Helmholtz coil)	housing	3 A/m	3 A/m	OK
Voltage dips, short-term interruptions	mains UT	<5% 0,5 per.	<5% 0,5 per.	OK
Test levels as % of UT duration in supply frequency periods	100 - 240V 50 / 60 Hz	40% 100 per.	40% 100 per.	OK
All testing at 100 and 240V each 50 Hz		70% 500 per.	70% 500 per.	OK
		<5% 5000 per.	<5% 5000 per.	OK

11 Disposal

 **Warning:** Ensure proper disposal of CardioSecur *ACTIVE* (electronic waste).

Labelling on CardioSecur *ACTIVE*, its packaging and this user manual indicate, that after its lifetime CardioSecur *ACTIVE* must not be disposed together with ordinary household refuse. Dispose of CardioSecur *ACTIVE* separately, so that harm to the environment or human health is not caused by uncontrolled disposal.

12 Labelling

www.cardiosecur.com
Personal MedSystems GmbH
Hansaallee 154
D-60320 Frankfurt am Main
Made in Germany

13 Contact

Personal MedSystems GmbH
Hansaallee 154
60320 Frankfurt am Main
Germany

Phone: +49-(0)69-907477-81
Fax: +49-(0)69-907477-84
E-Mail: info@cardiosecur.com
Website: www.cardiosecur.com

HRB 114144 B
USt.-Id.-Nr.: DE 260540641

 **Information:** Refer to our customer support in case of queries regarding CardioSecur *ACTIVE*
by phone: +49-(0) 69-907477-81
or via e-mail: info@cardiosecur.com

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